



CL 299 presentation 5

LNP Manufacturing and Process Simulation

Presented by:

Rhythem Soni 24110296

Aman Bola 24110032

Saket Jonwal 24110311

Supervisor:

Prof. Kaustubh S. Rane

Department of Chemical Engineering

Energy Balance

Equipment	Typical Lab Power Range
Syringe pump	5–20 W
Magnetic stirrer	5–15 W
Peristaltic pump	20–60 W
TFF circulation pump	80–120 W

$P_{\text{total}} = \mathbf{200\text{ W}}$

Energy consumption of the system is calculated using the equation:

$$E = P \times t$$

$$E = 0.22\text{ kWh}$$

What is in the Code?

- A Python-based simulator that takes feed compositions and flow rates as input and automatically calculates every stream in the process — including TFF splits and ethanol removal.

What You Input in the Code

- **Molar Ratios — Organic Phase**

Ionizable Lipid

Helper Lipid Cholesterol

PEG Lipid

Cholesterol

Ratios converted to wt% using MW of each component

- **Molar Ratios — Aqueous Phase**

mRNA

Buffer

Absolute concentrations set separately as mg/mL inputs

- **Feed Flow Rates**

Aqueous feed flow (mL/min)

Organic feed flow (mL/min)

Changing these two values changes every downstream flow automatically

- **Feed Concentrations**

mRNA conc in aqueous feed (mg/mL)

Buffer salt conc (mg/mL)

Total lipid conc in organic feed (mg/mL)

- **Process Parameters**

Encapsulation efficiency (fraction, e.g. 0.90)

Target ethanol vol% (e.g. 8.0%)

TFF Volume Reduction Ratio — VRR

Diavolumes N (for $C = C_0 \cdot e^{-N}$)

LNP retention at membrane

Salt sieving coefficient

Sterile filter dead volume (mL)

Process parameters are fixed equipment settings — not changed per run

Key Equations Used in the Model

Mole Fraction → Mass Fraction

$$wf_i = (x_i \times MW_i) / \Sigma(x_j \times MW_j)$$

Diafiltration Ethanol Removal ($C = C_0 \cdot e^{-N}$)

$$C_{\text{ethanol}} = C_0 \times \exp(-N)$$

Standard diafiltration washout equation. N = number of diavolumes. $N=5 \rightarrow 99.3\%$ removal. Used in FDA submissions and industry protocols.

Dilution Buffer Calculation

$$V_{\text{dil}} = V_{\text{EtOH}} / (\text{target\%} / 100) - V_{\text{S12}}$$

Back-calculates dilution buffer volume required to reduce ethanol from post-mixer level to target vol%. Adjusts automatically for any FRR.

Buffer Salt Removal with Sieving Coefficient

$$C_{\text{salt}} = C_0 \times \exp(-\sigma \times N)$$

Extension of ethanol model. σ (sieving coefficient) = fraction passing per diavolume. For 100 kDa membrane + small salts, $\sigma \approx 1.0$.

TFF Volume Reduction Ratio

$$VRR = V_{\text{S13}} / V_{\text{retentate}}$$

Sample Output (aq = 9 mL/min, org = 3 mL/min)

Stream	Flow (mL/min)	Ethanol (vol%)	LNP Mass Flow (mg/min)
Aqueous Feed	9.00	—	2.70 (mRNA)
Organic Feed	3.00	98.7%	30.00 (lipids)
Crude LNP (S12)	12.00	24.7%	32.70
Dilution Buffer	25.02	—	—
Diluted LNP (S13)	37.02	8.0%	32.70
DF Buffer → TFF	37.02	—	—
TFF Permeate	66.64	~99%→0	1.65 (lost)
TFF Retentate	7.40	0.2%	29.28
Final Product	7.40	0.2%	29.28

Balance Error
0%

Encapsulated mRNA Recovery
89.1%

Ethanol Removed (N=5 diavolumes)
99.33%

What the Simulator Outputs

- **Feed Streams**
Flow rate (mL/min)
Concentration of each component (mg/mL)
Mass flow of each component (mg/min)
- **Final Product Stream**
Final product flow rate (mL/min)
Concentration of all components in product
Mass flow of each component in final product
Encapsulated mRNA mass flow and recovery %
- **Calculated Process Flows**
Dilution buffer added (auto-calculated from target ethanol %)
Ethanol vol% at each stage
DF buffer flow (auto-scaled from diavolumes × retentate)
TFF permeate flow (auto-calculated from VRR)
- **Balance Verification**
Total mass IN vs total mass OUT
Volumetric flow accounting (IN = OUT)
Balance error (%)
Sterile filter dead volume noted separately

PAT

PAT classification - INLINE

Now we have classification of PAT techniques used in LNP manufacturing.

- Dynamic light scattering(DLS)
measure particle size and distribution & PDI (<0.2).
 - Zeta potential analysis
measure surface charge that indicates the stability of suspension.
 - Optical turbidity measurement(turbidity probes)
it detects incorrect lipid mixing
aggregation and wrong molecule formation.
 - PVM
capturing image while LNPs are forming
 - NIR
monitor composition and concentration of components in real time without sampling.
-

- UV–Vis Spectroscopy
it monitors Concentration of lipids, Encapsulated drug / RNA, Particle formation kinetics and Impurities.
- pH sensor
at the buffer mixing for correct pH.
- flow meter
measure flow rate of liquid streams.
- temp prob
- conductivity sensor
buffer dilution verification.
- Pressure sensor
in TFF filtration.
- Raman spectroscopy

Online

- High-Performance Liquid Chromatography (HPLC)
On-line possible via auto-samplers for drug loading, encapsulation efficiency, or degradation monitoring.
- Raman spectroscopy
used as a downstream analytical tool to measure how much mRNA is successfully encapsulated inside lipid nanoparticles.
- Charged aerosol detectors (CAD)
Used with chromatography to detect and quantify lipid components in LNP formulations that do not absorb UV light
- Refractive index (RI)
Measures changes in refractive index to detect lipids and polymers in LNP samples during chromatographic analysis.

- ultraviolet(UV) diode array
Detects UV-absorbing molecules such as RNA or drugs in LNP formulations while simultaneously recording their absorption spectra
- LC-MS
lipid identification and impurity detection

Atline

- Raman spectroscopy
- Zeta potential analysis
- DLS

- NMR
lipid structure and molecular interactions
- Cryo-TEM
Visualize morphology and internal structure of LNPs
confirm nanoparticle formation and observe structure of lipid core.

Resources:

- M. Mehta and T. A. Bui, “Lipid-Based Nanoparticles for Drug/Gene Delivery: An Overview of the Production Techniques and Difficulties Encountered in Their Industrial Development,” *ACS Materials Au*, vol. 3, no. 6, pp. 600–619, 2023, doi: 10.1021/acsmaterialsau.3c00032.
- J. Auclair and A. S. Rathore, “Analysis of Lipid Nanoparticles,” *LCGC North America*, vol. 41, no. 6, pp. 216–219, Jun. 2023, doi: 10.56530/lcgc.na.ii3372s8.
- S. Jin and E. Reverdatto, “Encapsulation of mRNA in Therapeutics Like Lipid Nanoparticles Probed by Deep-UV Resonance Raman Spectroscopy,” *Analytical Chemistry*, vol. 97, 2025, doi: 10.1021/acs.analchem.5c04246.
- J. A. Kulkarni, D. Witzigmann, S. Chen, P. R. Cullis, and R. van der Meel, “On the formation and morphology of lipid nanoparticles containing ionizable cationic lipids and siRNA,” *ACS Nano*, vol. 15, no. 5, pp. 9078–9093, May 2021, doi: 10.1021/acsnano.1c03058.
- https://drive.google.com/drive/folders/10AyL6bArZ2kjsxkyHl5krHhYheAy40hMO?usp=drive_link
- D. W. Green and M. Z. Southard, *Perry's Chemical Engineers' Handbook*, 9th ed. New York, NY, USA: McGraw-Hill Education, 2019.

- G. Towler and R. Sinnott, *Chemical Engineering Design: Principles, Practice and Economics of Plant and Process Design*, 3rd ed. Oxford, U.K.: Elsevier, 2022.
- N. M. Belliveau, J. Huft, P. J. Lin, S. Chen, A. K. K. Leung, J. B. Leaver, A. W. Wild, J. B. Lee, R. J. Taylor, Y. K. Tam, C. L. Hansen, and P. R. Cullis, “Microfluidic synthesis of highly potent limit-size lipid nanoparticles for in vivo delivery of siRNA,” *Molecular Therapy – Nucleic Acids*, vol. 1, no. 8, p. e37, Aug. 2012, doi: 10.1038/mtna.2012.28.
- R. van Reis and A. L. Zydney, “Membrane separations in biotechnology,” *Current Opinion in Biotechnology*, vol. 18, no. 3, pp. 208–211, Jun. 2007.

THANK YOU!!!